QUICK TUTORIALS

SPECTR-O-MATIC TOOLBOX

TABLE OF CONTENTS

Create specdata objects

Create specdata from variables

Create specdata arrays

Create specdata array from a matrix

Load data from a text file

Load multicolumn text file

Load multiple files

Multiple multicolumn files

Operate with spectra

Scalar operations

Operations with spectra and scalars

Arithmetic operations on spectra

Calculate sum and mean

Bin and smooth

Manipulate X axis

Shift and join spectra

Select and search

Select by direct indexing

Search for keywords

Find index

Automatic keyword index

Organize data by categories

Categorical arrays

Split-apply-combine workflow

Categorical index

Split-apply operations

Create parameter table

<u>Use external categorical index</u>

CREATE SPECDATA OBJECTS

SPECTR-O-MATIC TOOLBOX

CREATE SPECDATA FROM VARIABLES

1. Create some X and Y arrays:

```
x = 0:0.1:pi;
y = sin(x);
```

2. Then create a spectrum and name it "sinx":


```
s = specdata(x,y,'sinx');
```

3. Now plot the spectrum:

```
figure; plot(s)
```

X	1x32 d	ouble					
0	0.1	0.2	0.3	0.4	•••	3	3.1
у	1x32 d	ouble					
0	0.1	0.199	0.296	0.389	•••	0.141	0.042

s 1x1 specdata



CREATE SPECDATA ARRAYS


4. Create a new Y variable:

$$y = cos(x);$$


5. Create a second spectrum:

6. Plot both spectra at once:

```
figure; plot(s)
```



s 1x2 specdata



CREATE SPECDATA ARRAY FROM A MATRIX

1. Y is a matrix:


```
x = 0:0.1:pi;
y = [sin(x); cos(x)];
```

2. Create spectra:


```
s = specdata(x,y,{'sinx','cosx'});
```

3. Plot the spectra:

```
figure; plot(s)
```



s 1x2 specdata



LOAD DATA FROM A TEXT FILE

data1.txt

Wavelength	CD
350	13.428
350.5	13.247
351	13.041


•••

Load file

```
d = specdata.load('data1.txt');
```

Plot

figure; plot(d)



LOAD MULTICOLUMN TEXT FILE

data2.txt

Wavelength	CD	HT	Abs
350	0.236	0.636	0.118
350.5	0.228	0.636	0.117
351	0.221	0.635	0.116

...

Load file

```
d = specdata.load('data2.txt');
```

Show properties ID, XType, YType

d 1x3 specdata

ID	XType	YType
		
'data2'	'Wavelength'	'CD'
'data2'	'Wavelength'	'HT'
'data2'	'Wavelength'	'Abs'

LOAD MULTIPLE FILES

data1.txt

Wavelength CD 350 -7.226 350.5 -7.139 351 -7.097

Load files

data2.txt

Wavelength	CD
350	-8.592
350.5	-8.579
351	-8.533

data3.txt


Wavelength	CD
350	3.113
350.5	3.061
351	3.046

```
d = specdata.load('*.txt');
```

Plot

figure; plot(d)

d 3x1 specdata



MULTIPLE MULTICOLUMN FILES

data1.txt

data2.txt

data3.txt

Wavelen	CD	Abs
350	-7.226	0.229
350.5	-7.139	0.227
351	-7.097	0.224

Wavelen	CD	Abs
350	-8.592	0.556
350.5	-8.579	0.552
351	-8.533	0.549

Wavelen	CD	Abs
350	3.113	0.500
350.5	3.061	0.498
351	3.046	0.495


.

Load files

Plot 2nd column of *d*

figure; plot(d(:,2))

d 3x2 specdata



OPERATE WITH SPECTRA

SPECTR-O-MATIC TOOLBOX


SCALAR OPERATIONS

Load data

load data.mat

```
>> whos
  Name Size Bytes Class
  a 3x1 39090 specdata
```

>> plot(a)



Extract Y values (at position X)

$$m = a.Yx(675);$$

>> m'

0.6360

0.5500

0.8230

Maximal values of Y (within a range)

$$mx = a.max([400 500])$$

>> mx'

1.1480

0.9040

1.1640

Find peaks (of a minimum amplitude)

$$p = a(1).peaks(0.1)$$

>> p


438.0000 1.1480

674.5000 0.6360

OPERATIONS WITH SPECTRA AND SCALARS

Divide spectra by scalar array

$$b = a / m;$$




Shortcut

$$b = a.norm(675);$$

Add scalar array to spectra


$$b = a + [0, 0.5, 1];$$



ARITHMETIC OPERATIONS ON SPECTRA


Subtract one spectrum from an array

$$b = a(2:3) - a(1);$$



Binary operations with arrays of spectra


$$b = -0.5*a + a^2;$$




CALCULATE SUM AND MEAN

Sum spectra

$$b = sum(a);$$



Average spectra




BIN AND SMOOTH

Bin data points (boxcar average)

b = a.bin(20);


>> [b.dim] 40 40 40



Smooth (moving average)

b = a.smooth(20);

>> [b.dim] 801 801 801




MANIPULATE X AXIS

Trim spectra (set X limits)

b = a.setxlim([600 700]);


>> [b.dim] 201 201 201



Change X axis (interpolate spectra)

b = a.setx(600:0.1:700);


>> [b.dim] 1001 1001 1001



SHIFT AND JOIN SPECTRA


Shift spectra along X axis

$$b = a.shiftx(-675);$$



Join spectra

```
b = a(1).setxlim([400 600]);
c = a(2).setxlim([500 700]);
d = merge([b,c]) + 0.5;
>> plot(b,c,d)
```



SELECT AND SEARCH

SPECTR-O-MATIC TOOLBOX

SELECT BY DIRECT INDEXING

Load data

```
load data.mat
>> whos
  Name Size Bytes
                     Class
        3x3
                     specdata
  dat
              117506
>> dat.get('ID')
   'SampleA fh'
                  'SampleA fv'
                                'SampleA iso'
   'SampleB fh'
               'SampleB fv'
                                'SampleB iso'
   'blank fh'
                'blank fv'
                                 'blank iso'
```

Select a single spectrum

A = dat(1);

```
>> A.ID
  SampleA fh
Select multiple spectra
A = dat([1,4,7])
>> A.get('ID')
                                   'SampleA iso'
                   'SampleA fv'
```

Select a row of spectra

'SampleA fh'

```
A = dat(3,:)
>> A.get('ID')
              'blank fv'
                             'blank iso'
   'blank fh'
```

SEARCH FOR KEYWORDS

Find spectra by ID

```
A = dat.find('SampleA');

>> A.get('ID')
    'SampleA fh'
    'SampleA fv'
    'SampleA iso'
```

Search for multiple keywords (match ANY)

```
A = dat.find({'SampleA', 'SampleB'});
>> A.get('ID')
   'SampleA fh'
   'SampleB fh'
   'SampleA fv'
   'SampleB fv'
   'SampleA iso'
   'SampleB iso'
```

Search different properties

Search for multiple keywords (match ALL)

```
A = dat.find('ID', 'SampleA', 'ID', 'iso');
>> A.get('ID')
    'SampleA iso'
```

FIND INDEX

Find index of spectra

```
a = dat.findindex('SampleA');

>> whos a
   Name Size Bytes Class
   a 9x1 9 logical
>> a'
   1 0 0 1 0 0 1 0 0
```

Select spectra using index

```
A = dat(a);
>> A.get('ID')
    'SampleA fh'
    'SampleA fv'
    'SampleA iso'
```

Combine results (Boolean search)

```
a = dat.fi('SampleA');
b = dat.fi('SampleB');
ab = a | b; % SampleA + SampleB
>> a'
                          0
                                         0
>> b'
           0
                0
                          0
                               0
       1
  0
                                         0
>> ab'
                          0
                                         0
i = dat.fi('iso');
```

```
i = dat.fi('iso');
abi = ab & i; % (SampleA + SampleB) * iso

>> i'
    0     0     0     0     0     1     1     1
>> abi'
    0     0     0     0     0     1     1     0 >
```

AUTOMATIC KEYWORD INDEX

Create an automatic keyword index

'SampleA iso'

```
x = dat.autoindex;

>> x
    SampleA: [9x1 logical]
    SampleB: [9x1 logical]
    blank: [9x1 logical]
        fh: [9x1 logical]
        fv: [9x1 logical]
        iso: [9x1 logical]

Select spectra using index

A = dat(x.SampleA);

>> A.get('ID')
    'SampleA fh'
    'SampleA fv'
```

Select spectra (Boolean OR)

```
A = dat(x.SampleA | x.SampleB);
>> A.get('ID')
    'SampleA fh'
    'SampleB fh'
    'SampleB fv'
    'SampleB iso'
```

Select spectra (Boolean AND)

```
A = dat(x.SampleA & x.iso);
>> A.get('ID')
    'SampleA iso'
```

ORGANIZE DATA BY CATEGORIES

SPECTR-O-MATIC TOOLBOX

CATEGORICAL ARRAYS

Load data

```
load data.mat
>> whos
  Name Size Bytes
                      Class
        3x3
  dat
              117506
                      specdata
>> dat.get('ID')
   'SampleA fh'
                  'SampleA fv'
                                 'SampleA iso'
   'SampleB fh'
                'SampleB fv'
                                 'SampleB iso'
   'blank fh'
                 'blank fv'
                                  'blank iso'
```

Create a categorical array (MATLAB 2016b+)

```
c = dat.catfind({'SampleA','SampleB'});
>> c
    9x1 categorical array
```

Select spectra using categorical array

```
A = dat(c=='SampleA');
>> A.get('ID')
    'SampleA fh'
    'SampleA fv'
    'SampleA iso'
```

SPLIT-APPLY-COMBINE WORKFLOW

1. Create categorical array of samples

```
s = dat.catfind({'SampleA', 'SampleB'});
```

2. Average spectra per group

```
g = findgroups(s);
A = splitapply(@mean, dat, g);
>> A
    2×1 specdata array
% splitapply ignores <undefined> samples
```

3. Find baseline spectra

```
B = dat.find('blank');
```

4. Split data and subtract baselines per group

```
f = @(x) \{x-B\}; % function <math>f(x) = \{x-B\}
A = splitapply(f, dat, g);
```

```
>> A
    2×1 cell array
    [3×1 specdata]
    [3×1 specdata]
```

5. Concatenate result (if needed)

```
A = cat(1, A{:});
```

```
>> A
6×1 specdata array
```

CATEGORICAL INDEX

Define variables

```
vars = struct;
vars.Sample = {'SampleA','SampleB'};
vars.Code = {'iso','fh','fv'};
```

Create index of variables

```
C = dat.catindex(vars);
>> C
                     Sample
                                 Code
                                 fh
   SampleA fh
                  SampleA
   SampleA fv
                  SampleA
                                 fv
   SampleA iso
                  SampleA
                                 iso
   SampleB fh
                  SampleB
                                 fh
   SampleB fv
                                 fv
                 SampleB
   SampleB iso
                  SampleB
                                 iso
   blank fh
                  <undefined>
                                 fh
   blank fv
            <undefined>
                                 fv
                <undefined>
   blank iso
                                 iso
```

Select spectra

```
A = dat(C.Sample=='SampleA');
>> A.get('ID')
    'SampleA fh'
   'SampleA fv'
    'SampleA iso'
A = dat(isundefined(C.Sample));
>> A.get('ID')
   'blank fh'
   'blank fv'
   'blank iso'
A = dat(C.Sample=='SampleA' ...
         & C.Code=='iso');
>> A.get('ID')
                                                  >
   'SampleA iso'
```

SPLIT-APPLY OPERATIONS

Average spectra grouped by 'Sample'


```
A = dat.splitop(@mean, C, 'Sample');
>> A
    3×1 specdata array
% splitop treats <undefined> as a separate group
```

Subtract baselines grouping by 'Code'

```
b = dat.fi('blank');
A = dat.splitbinop(@minus,~b,b,C,'Code');
>> A
    1×6 specdata array
```

Plot spectra by sample group

```
C = A.catindex(vars);
A.splitop(@plotf, C, 'Sample');
```



CREATE PARAMETER TABLE

1. Create categorical index

>> C

		<u>Sample</u>	<u>Code</u>
SampleA	fh	SampleA	fh
${\tt SampleA}$	fv	SampleA	fv
${\tt SampleA}$	iso	SampleA	iso
SampleB	fh	SampleB	fh
SampleB	fv	SampleB	fv
SampleB	iso	SampleB	iso

2. Add a calculated parameter column

$$C.Amax = A.max([600 700]);$$

>> C

	<u>Sample</u>	<u>Code</u>	Amax
SampleA fh	SampleA	fh	0.52855
SampleA fv	SampleA	fv	0.51194

3. Pivot table with 'Code' in columns

P = unstack(C, 'Amax', 'Code');

>> P
 <u>Sample</u> <u>iso</u> <u>fh</u> <u>fv</u>
 SampleA 1.0947 0.5286 0.5119
 SampleB 0.8120 0.3831 0.4333

4. Pivot table with 'Sample' in columns

P = unstack(C, 'Amax', 'Sample');

>> P

<u>Code</u> <u>SampleA</u> <u>SampleB</u>

fh 0.52855 0.38317

fv 0.51194 0.43329

iso 1.0947 0.81202

USE EXTERNAL CATEGORICAL INDEX

DataIndex.xls

FileName	Sample	Code
data3	SampleA	iso
data1	SampleA	fh
data2	SampleA	fv

. . .

Load data files

```
dat = specdata.load('data*.txt');
>> dat.get('ID')
   'data1'
   'data2'
   'data3'
   ...
```

Synchronize external index table to data

Synchronize preloaded index table to data